

Overnight Oximetry – Sleep Study

Step 1:

At your normal bedtime, put the probe on a finger with no rings or nail polish. Tape the probe in place to keep it there as shown in picture. DO NOT TURN OFF UNTIL TEST COMPLETE IN MORNING.




Step 2:

Turn the machine on by pressing the  on button. When there is a good signal, a green light flashes with each heartbeat.

If you wake up during the night to go to the bathroom or probe feels too tight, you can put it on a different finger with no rings or nail polish.

Step 3:

In the morning turn the machine off by pressing and holding the  on button until the machine turns off.

Step 4:

Bring the machine back to Transcend Medical or Return by Mail with the enclosed envelope. Sign Enclosed AOB where highlighted and include in package. (Thanks)

Please let us know if you have any problems getting the machine back in a timely manner.

You can call us at 800-403-3740

www.transcendmedical.net

**133 Woods Cove Road
Scottsboro, AL 35768
35976**

**2001 Henry St
Guntersville, AL**