

# OxyGO® FAQs

## FREQUENTLY ASKED QUESTIONS

### IS AN OXYGO THE SAME AS AN INOGEN UNIT?

Yes, OxyGo portable oxygen concentrators are the same devices as Inogen portable oxygen concentrators. The only difference is the brand name and that you will be supported by your local provider.

OxyGo NEXT®	=	Inogen One G5®
OxyGo FIT®	=	Inogen One G4®

### WHERE ARE OXYGO CONCENTRATORS MADE?

OxyGo POC's are made in the same manufacturing facility, using the same components, as Inogen portable oxygen concentrators.

### WHY SHOULD I CHOOSE OXYGO?

OxyGo POCs are provided and supported by local providers in your area. Your local provider will be able to deliver fast local service for your concentrator while giving you the one-on-one treatment you require.

### WHAT COMES WITH OXYGO?

Each OxyGo comes with everything you need to Keep Going! Each unit comes with a single battery, car charger, AC/DC wall charger, cannula, quick start guide and carrying case. There are also optional upgrade packages with extra batteries and double batteries available.

### CAN I FLY WITH OXYGO?

Yes! All OxyGo POCs comply with the applicable rules and regulations for domestic and international air travel. The FAA has established acceptance criteria for POCs used on aircraft. Each unit has markings on the back label indicating that they are compliant. We're happy to announce that our OxyGo POCs are compliant, and bear the necessary labeling on the back of the device.\*

\*It is always recommended to contact your individual airline to verify their requirements when travelling with a POC.

BATTERY LIFE\*

OxyGo NEXT™

OxyGo FIT™

Setting	Output (ml/min)	Single Battery (hr/min)	Double Battery (hr/min)	Single Battery (hr/min)	Double Battery (hr/min)
1	210	6:35	13:28	2:38	4:54
2	420	5:03	10:22	2:14	4:29
3	630	3:27	7:01	1:23	2:40
4	840	2:25	4:55		
5	1,050	1:51	3:48		
6	1,260	1:26	2:58		

\*Battery life statistics based on average values may vary up to +/- 10% depending on patient use.



Keep Going.™