

Flow Rate LPM	Mini M-6 Cylinder	ML-6 Cylinder	“C” Cylinder	“D” Cylinder	“E” Cylinder
.5	16.5 hours	17.2 hours	24.2 hours	41.9 hours	68.9 hours
.75	11.0 hours	11.4 hours	16.2 hours	28.0 hours	45.9 hours
1.0	8.3 hours	8.6 hours	12.1 hours	21.0 hours	34.4 hours
1.5	5.5 hours	5.7 hours	8.1 hours	14.0 hours	23.0 hours
2.0	4.1 hours	4.3 hours	6.1 hours	10.5 hours	17.2 hours
2.5	3.3 hours	3.4 hours	4.9 hours	8.4 hours	13.8 hours
3.0	2.8 hours	2.9 hours	4.0 hours	7.0 hours	11.5 hours
3.5	2.4 hours	2.4 hours	3.5 hours	6.0 hours	9.8 hours
4.0	2.1 hours	2.1 hours	3.0 hours	5.2 hours	8.6 hours
5.0	1.7 hours	1.7 hours	2.4 hours	4.2 hours	6.9 hours
6.0	1.4 hours	1.4 hours	2.0 hours	3.5 hours	5.8 hours

Note: Usage times for conserving devices vary depending upon cylinder size, device type, and patient. This chart is meant to be a general guide only; your actual usage times may vary.